



ibis WHEELCHAIR TO CAR ACCESS

The IBIS® Wheelchair to Car Access has been designed and manufactured in Australia by Para Mobility. This is a transformative device that allows a person with a disability (PWD) and mobility limitations the option of getting in a car, SUV, 4WD or van.

The IBIS is designed to attach to the hinges of the vehicle door fully eliminating the need for any permanent modifications and giving users the freedom to travel in multiple vehicles. Our system is conveniently powered by the car's 12 volt auxiliary outlet and can be operated by the carer via the hand control. With a simple button press the carer gently guides the participant into the vehicle without any manual lifting.

By simplifying vehicle access, the IBIS allows users to travel more freely, whether for daily errands, social outings, or medical appointments. This increased mobility enhances overall quality of life.

FEATURES

- Fully portable and compatible with most vehicle makes and models
- No vehicle modification required
- Operation is by the carer with no manual lifting required
- Press button hand control powered by the vehicles 12 volt auxiliary outlet (optional battery pack available)
- Easily removed and stored in the boot when travelling
- Lightweight for easy handling – 17kgs (assembled in 4 parts)
- Retains its resale value, making it a smart long-term investment
- Two (2) size options suitable for small/medium and large vehicles
- TGA approved & NDIS confirmed

INCLUDES

- Bracket
- Lifting Frame (Frame, Extension Bar and Chest Piece)
- Hand Controller with 12V power cable
- Sling
- Carry/Storage Bag



Freedom to explore



Para Mobility
Making everyday living easier

Phone 1300 444 600 +61 2 9651 4446
Delivery, installation + servicing Australia wide
enquiries@paramobility.com.au
www.paramobility.com.au

LEARN MORE

