



ParaMobility

Making everyday living easier

Phone 1300 444 600

Unit 31/276 New Line Rd

Dural NSW 2158 Australia

enquiries@paramobility.com.au

www.paramobility.com.au


Pelican Pool & Spa Hoist Measuring Guide

Customer Name _____ Phone _____

Pool Name _____

Install Address _____

Pelican Rim Setback (Aerial View)



Pelican setback is the standard setback distance (100mm).

IMPORTANT NOTICE: WATER DEPTH

In accordance with the Disabilities Standards Code, Disability (Access to Premises – Buildings) Standards 2010 Part D5, the Pelican base rim must be positioned on the pool or spa surrounds where the **water depth is not less than 500mm and not more than 1300mm.**

Pool depth at the point of entry _____ mm Pool depth at location _____ mm

Which profile best matches your pool? ☐ A ☐ B ☐ C ☐ D ☐ E

Please include measurements in the profile you selected (see page 2).

Do you have the minimum concrete width & thickness to mount the hoist on pool deck? ☐ Y ☐ N

Which Dulux colour would you like?



☐ Pearl White Gloss

☐ Charcoal Metallic Pearl Gloss

☐ Silver Kinetic Pearl Satin

☐ Safety Yellow Gloss

☐ French Blue Gloss

☐ Primrose Gloss

What is the base substrate material the hoist will be installed on? _____
(Tile, concrete, timber, etc.)

Property access considerations/risks? _____



ParaMobility

Making everyday living easier

Phone 1300 444 600

Unit 31/276 New Line Rd

Dural NSW 2158 Australia


enquiries@paramobility.com.au

www.paramobility.com.au

Pelican Pool & Spa Hoist Measuring Guide

ROUND MEASUREMENTS TO THE NEAREST 10MM

Profile A




A1 COPING WIDTH

A2 INTERNAL RISER

A3 EXTERNAL RISER


NOTE: Measurement 2 may be different to measurement 3.

Profile B



B2 INTERNAL RISER


Profile C



C1 GRATE WIDTH


C2 INTERNAL RISER

Profile D



D1 GRATE WIDTH

Profile E



E1 COPING WIDTH

E2 INTERNAL RISER

E3 EXTERNAL RISER

E4 WIDTH OF
SPA SEAT

NOTE: Measurement 2 may be different to measurement 3.